

# Aquilos2 FIB User Guidelines ABC-061

## Liu Box – Aquilos FIB

**Main Contact:** Jenny [jennifer.aronson@yale.edu](mailto:jennifer.aronson@yale.edu) 203-500-5012  
**Training/Technical:** Seven [wangbiao.guo@yale.edu](mailto:wangbiao.guo@yale.edu) 203-390-0472

chris.olsen@thermofisher.com + john.petrus@thermofisher.com + Yale Vendor ID, Trainer 12/2021 ken.wu@thermofisher.com

### 1. Non-Liu Lab users: [User Request Email Form](#)

Jenny: emails & access: ID, Box-Aquilos FIB.

### 2. All users: [Usage Log](#)



### 3. Liquid Nitrogen

- Notify Jenny **before** tank empty.
- To switch LN2 tanks:
  1. Close 'Vent' and 'Liquid' valves on empty tank.
  2. Move hoses from empty tank to new tank, connect to correct outlets (Vent - Vent, Liquid - Liquid).
  3. Open valves on new tank fully.
  4. Move empty tank to room 059C, label with paper marked 'RETURN'.

## Instructions: Liu Box – Aquilos FIB – Guidelines, Access, Trainings

### 1. Before operation

- 1) Check preparation pot at **purge** state, turn on **heat** button to make sure tools dried
- 2) Check Chamber Pressure **< 4e-4Pa**, write exact value in [Usage Log](#)
- 3) Check **Beam shift** is zero
- 4) If nobody used Aquilos2 one day prior, perform sputtering and clean needle for 2 minutes (skip if not needed).
- 5) Home stage.



### 2. Sputtering (1<sup>st</sup> time, optional. **No sample loaded at this point.**)


- 1) Prepare for sputtering (**30 mA, 10 Pa, 15 sec**)
- 2) Run
- 3) When stage stops moving, click **Recover from sputtering**.

### 3. Clean needle (optional)

- 1) Change Flow Duration to 2 mins
- 2) **Purge**
- 3) For this time, no need to click Grid1 or Grid2.


### 4. Cool down

- 1) After preparation steps, start to cool down stage and shield
- 2) Turn on Flow DDE software
- 3) Communication → Open communication
- 4) Turn on FlowView software
- 5) Adjust flow rate to 200 mg/s
- 6) Wait until cryo-stage temperature < -190 °C.




## 5. Load samples

- 1) If unsure, ask Seven for training
- 2) After loading sample, turn on E-beam, choose Grid1 or Grid2
- 3) Go to Mapping position
- 4) Check sample quality, such as ice, thickness, etc.
- 5) Decide if experiment will continue or which grid to use.
- 6) If good, link Z position (While E beam is on) Do this for both grid 1 and 2
- 7) Change detector to ETD, Mode 1: standard



## 6. Montage

- 1) Turn on software Maps 3.18
- 2) Create project name, choose project path... Save snapshots, Maps, AutoTEM data at: **/ShareData/User\_Snapshots/year-m-d\_User\_SampleName(Collaborator)** (e.g. /ShareData/User\_Snapshots/2022-01-03\_Seven\_Algae)
- 3) Change magnification on electron beam ~100x, take snapshot image on Maps. If two grids, move to other grid and take snapshot again.
- 4) Add tile set (choose 'seven' template):  
Set tile parameters: Tile HFW: 550  $\mu\text{m}$ , Resolution: 3072\*2048, Dwell time: 1  $\mu\text{s}$ , Frames: 1, Focus Strategy: None, Contrast Brightness: First tile, Focus: First tile, Order: Spiral Out
- 5) Click Run
- 6) When finished, save snapshot image to local
- 7) Adjust eucentric positions manually (rather than AutoTem focus) to check Z-Y link, align ion beam and E beam positions when too far apart. Click ctrl + F12, then F6, left click to drag to desired position. Add lamella position (right side, 2<sup>nd</sup> panel of UI) after storing angle.



After aligning all lamella positions, **uncheck** the Z-Y link.

- 8) To use FIB data or Cryo-CLEM images on SPC, login in browser and save to your Yale Box account (or temporary: Liu [Box](#) Aquilos FIB – FIB\_images. Move/delete your data from this folder after FIB session).
- 9) Save all data (snapshots, Maps\_Project, AutoTEM\_Project) to: ShareData/User\_Snapshots – Data\_User\_Sample (e.g. 2022-01-03\_Seven\_Algae). Folder saved monthly to our clusters. Never use USB drive on SPC (see note: No USB Allowed).

[alternate: Seven Box Gmail: Aquilos2\_User\_Temp\_Data: <https://app.box.com/s/3jv1pxvdqbg1apmqwoxtjc5tn76db78x>]

## 7. Import CLEM image (if needed)

- 1) File → Import images
- 2) Choose Tile Set → Alignment
- 3) Two options: Choose Align or Fine Alignment.

## 8. Choose lamellae

- 1) Choose good and clean areas
- 2) Right click, select set lamella site here.

## 9. Sputtering (2<sup>nd</sup> time, optional)

- 1) Prepare for sputtering (**30 mA, 10 Pa, 15 sec, No1**)
- 2) Run
- 3) When stage stops moving, click Recover from sputtering.

## 10. GIS coating

- 1) Select Grid1 or Grid2
- 2) change Flow Duration, this time depends on samples.
- 3) Remember this time and modify based on your milling results.

## 11. Sputtering (3<sup>rd</sup> time, optional)

- 4) Prepare for sputtering (**30 mA, 10 Pa, 15 sec**)
- 5) Run
- 6) When stage stops moving, click Recover from sputtering.

## 12. AutoTEM Cryo2.2

- 1) Turn on software AutoTEM Cryo2.2
- 2) Make sure Maps software is on, choose project "on-line", then AutoTEM and Maps link. To choose more lamellae on Maps, AutoTEM can read these simultaneously.
- 3) Choose template → Apply or Apply to all
- 4) Choose each lamella position. Click "go to" in the UI, then click "Update"
- 5) At this point, it is not necessary to adjust the stage position even if positions are not quite right. When all sites are finished, click "Run". Each site takes ~3 seconds.
- 6) Adjust the stage positions of E beam and ion beam:
  - a) At UI, go to Ebeam, make sure target is centered. Image is clear (adjust focus and stigmatism, WD 7 mm helps).
  - b) Adjust ion beam position to align with E beam. When finished, adjust magnification to 10,000 X. Adjust the milling angle to a larger number (milling angle 16 degrees...change back to 8 or 9 later). Return to AutoTEM, click "Update".
- 7) At autoTEM, adjust rectangles to milling positions.
- 8) Repeat for all sites.
- 9) May need to click the square ("stop") in UI to manually move to next position.
- 10) At this point, autoTEM only suggested for rough milling, especially for cell samples, two templates:
  - a) Seven bacteria fine milling
  - b) Seven cell rough milling

## 13. Preparation (AutoTEM, not reliable)

Helps to align eucentric position:

- 1) Maximal tilt step: default 15°, can lower to 10° if sample too thin
- 2) Milling angle: default 10° (lamella surface angle) v. stage angle + 7°
- 3) Choose Automatic (with manual fallbacks)
- 4) When finished, check ion beam and electron beam to align
- 5) **If beam shift not huge, use shift + left drag. If beam shift huge, double click to center area of interest in E beam, then make sure ✓ Z-Y link is checked, use stage Z to move (ctrl+F12) to center the area of interest on ion beam also. Smaller ion beam current recommended in case sample damaged during adjustment.**
- 6) Default WD on ion beam 19.0 mm, electron beam 7.0 mm. Can manually change back to help adjust focus
- 7) When finished, manually adjust to focus well at lamella site, ideally in middle of image.
- 8) **If position and focus are adjusted, click "update" in AutoTEM cryo software.**
- 9) On AutoTEM Cryo software, adjust lamella position, then click "Continue"
- 10) System will move to next lamella automatically.

## 14. Rough milling

- 1) After preparation, start rough milling. Milling current is sample to sample, best to modify parameters carefully.
- 2) At STRESS RELIEF CUTS, Milling Current: 0.5 nA, DCM rescan interval: 120s – 300s
- 3) At ROUGH MILLING, Milling Current: 0.5 nA, Pattern Type: Rectangle, DCM rescan interval: 120s – 300s
- 4) At MEDIUM MILLING, Milling Current: 0.5-0.3 nA
- 5) At FINE MILLING, Milling Current: 0.5-0.3 nA
- 6) At FINER MILLING, Milling Current: 0.1 nA
- 7) Click Finer Milling-Electron Image to take snapshot after rough milling.

E beam detector change to T2, mode 3: optitilt

## 15. Fine milling

- After AutoTem is finished, click to go to milling position.
- **Change angle to a smaller angle (e.g. 8 or 9).**
- Change T2 **and optitilt** (Mode 3)
- Polish lamella manually
- Suggested milling current 0.5 nA.
- Resolution (ion beam -> 3042 x 2048)
- 300 ns

- E beam (resolution 3042x 2048)
- 3  $\mu$ s
- Stop at ~600 nm.
- Use 50 pA to polish lamella to < 200nm.
- Do not prolong experiments, finish before 8:00 p.m. (vacuum system not ideal).

#### 16. Save samples

- 1) **Make sure the sample preparation pot has cooled down 30 minutes before the last lamella thinning.**
- 2) Take sample out of chamber
- 3) Turn N2 flow rate to 20 mg/s
- 4) Turn beam shift to zero
- 5) Turn off ion beam and electron beam
- 6) **Usage Log:** Record details in all fields and log out of Box/websites used on SPC.
- 7) **Grid Storage System Log:** Save milled FIB sample in puck 11 with box color, box name, sample name, session.
- 8) Use "sleep source" if no one using FIB 2+ days.

**Do not edit below this line...**

---

Seven 12/21:

Save FIB images to **server**: FIB\_images/[folder with your name] / \_\_\_\_\_

Save and check lamellae in Box Sync software on FIB computer (right): Computer/Data(D:)/ShareData/BoxSync/FIB\_images.

[Zoom training recordings 12/21

<https://yale.box.com/s/ngtbos2sskas4ui5r5bftnyhp5nhxwss>

<https://yale.box.com/s/u9oqu7ao0qrr7lkf3wkhpf42lh266d1>

<https://yale.box.com/s/e7oiasnnl7iq7khj5btuvk3u2yespw17>

<https://yale.box.com/s/tocwyxjyxr37ri6lqn77jttgcijsmurt>

<https://yale.box.com/s/8l8qlcasx5nzos71vu8svui79rnvhqx3>

Meng 11/21:

- All FIB SEM images on computer in FIB room: /desktop/shareddata/FIB images
  - ✓ All FIB experiment records, equipment operations, expense in Box.: Liu\_CryoET\_Lab/Aquilos\_FIB.
  - Training videos in Box and: @Gunpowder: /data/disk3/FIB/FIBTrainingVideo.
-